



PROTOS[®] LATEX

Examination Gloves



The **BEST** Gloves
For our **HEALTH**

SCAN ME



www.protosindonesia.com @protosindonesia



Lightly Powdered (AKD 10903410022)
Powder-Free (AKD 10903410021)

QUALITY STANDARD

- Marketing License for Medical Devices from The Minister of Health
- Conforms to ASTM D3578-05
- Manufactured under the ISO 9001:2008, Quality Management System

GLOVES SIZE

- Extra Small
- Small
- Medium
- Large
- Extra Large

FUNCTIONAL BENEFIT

- Protection from unwanted or dangerous substances
- Easy donning and helps prevent roll back
- Softness provides superior comfort and natural fit
- Beaded cuff makes donning easy
- Ambidextrous and straight fingers
- Eliminate powder induced irritation and dermatitis to the users*
- Textured surface on palm and fingers provide a better grip*
- Low level of extractable proteins and chemical residue, hence lower incidences of allergic reaction*

PRODUCT SPECIFICATIONS

TYPE Lightly Powdered, Powder-Free (Non-Sterile)
MATERIAL Natural Rubber Latex

DESIGN & FEATURES

	DISPOSABLE	AMBIDEXTROUS	BEADED CUFF	SURFACE	COLOR	POWDER
Lightly Powdered	✓	✓	✓	Smooth	Natural White	Absorbable Dusting Powder
Powder Free	✓	✓	✓	Micro Textured	O-White to Yellow	✗

STORAGE in dry condition at temperature not higher than 30° C
SHELL LIFE 5 years from the date of manufacturing
PACKAGING 10 boxes/case; 100 pieces/box (by weight)



PT. MENARA MEDIKA PRATAMA

Jl. RC Veteran No.9,
Kel. Bintaro, Kec. Pesanggrahan
Jakarta Selatan - 12330

Telp +62 - 21 - 29511968
Fax +62 - 21 - 73889661
Email marketing@protosindonesia.com