



PROTOS[®] STERILE LATEX

Surgical Gloves



The **BEST** Gloves
For our **HEALTH**

SCAN ME



www.protosindonesia.com @protosindonesia



Powdered (AKD 11603410258)
Powder-Free (AKD 11603410183)

QUALITY STANDARD

- Marketing License for Medical Devices from The Minister of Health
- Conforms to ASTM D3578-09e1
- Manufactured under the ISO 9001:2008, Quality Management System
- Sterilized by EO Gas

GLOVES SIZE

- 6.0
- 6.5
- 7.0
- 7.5
- 8.0

FUNCTIONAL BENEFIT

- Extra strength provides additional protection from surgical debris
- Easy donning and helps prevent roll back
- Softness provides superior comfort and natural fit
- Beaded cuff makes donning easy
- Fully anatomical design to reduce hand fatigue
- High strength and elasticity
- Micro textured surface provide excellent wet and dry grip
- Low level of extractable proteins and chemical residue, hence lower incidences of allergic reaction*

PRODUCT SPECIFICATIONS

TYPE Powdered, Powder-Free (Sterile)
MATERIAL Natural Rubber Latex

DESIGN & FEATURES

	DISPOSABLE	ANATOMIC	BEADED CUFF	SURFACE	COLOR	POWDER
Powdered	✓	✓	✓	Micro Textured	O-White to Yellow	Absorbable Dusting Powder
Powder Free	✓	✓	✓	Micro Textured	O-White to Yellow	✗

STORAGE SHELL LIFE Store in dry condition at temperature not higher than 30°
5 years from the date of manufacturing



PT. MENARA MEDIKA PRATAMA

Jl. RC Veteran No.9,
Kel. Bintaro, Kec. Pesanggrahan
Jakarta Selatan - 12330

Telp +62 - 21 - 29511968
Fax +62 - 21 - 73889661
Email marketing@protosindonesia.com